


**Leg Bags
Leg Bag Accessories
Drainage Bags**



More Choice, more comfort



Attention to Detail. Attention to Life.

Leg Bags, Leg Bag Accessories & Drainage Bags

More choice, more comfort

Caring for the details that make all the difference.

As a healthcare professional, you develop a bond of trust with those you help and their families – throughout their lives.

The role of Hollister is to support you in these vital relationships by focusing on the details. Experience has taught us that sometimes the smallest things make the greatest difference. You'll find this approach reflected in the design of all our products, our services and everything we do.

For more than 80 years, Hollister has been making a difference where it matters, playing a role in making life more rewarding and dignified. As we continue to develop products and services, our goal remains unchanged – to promote independence and quality of life.



Leg Bags

Hollister Leg Bags offer flexibility to suit a wide range of patient needs.

This discreet collection and removal solution complements an active lifestyle and blends user-friendly features with the optimum standards of hygiene.

- 500ml, 540ml and 800ml capacities
- Tubing lengths: 10cm or 50cm
- Waterproof material limits unwanted odours
- Comfort backing
- Anti-reflux valve helps protect against leakage and urinary tract infections caused by urine backflow
- Outlet valve features easy-to-handle lever tap
- Protected sampling port helps minimise needle stick injuries



'Christmas tree' adapter is designed to fit smoothly, while providing maximum security



Easy-to-handle twist-action lever tap

Leg Bag Straps

- Comfortable and cut to fit
- Anti-slip silicone bands
- Attachment to thigh or calf

Deluxe Leg Bag Straps

- Four different lengths
- Wider and more comfortable
- Moisture-absorbent
- Easy to use

Leg Bag Holsters

As an alternative to Leg Bag Straps, Leg Bag Holsters offer your patient the freedom to choose the most secure, discreet and comfortable option to suit their lifestyle.

- Available in three sizes each, for both thigh and calf
- Allow freedom of movement
- Colour coded for easy size identification

Night Drainage Bags

- 120cm inlet tubing – allows freedom of movement
- 2000ml capacity
- Anti-reflux valve protects against leakage and infection caused by urine backflow
- Drainable or non-drainable formats
- Drainable version features easy-to-handle lever tap and sampling port



Leg Bag Straps




Deluxe Leg Bag Straps



Leg Bag Holster




Night Drainable Bags
Drainable (sterile)




Night Drainable Bags
Non-drainable (non-sterile)


There's more to Hollister Continence Care than ever before

More choices. More solutions. More from the heart of Hollister.

Urinary Leg Bags		Order No.	Pack Size
	Sterile Leg Bag 540ml without extension tube	9814	10
	Sterile Leg Bag 500ml 10cm inlet (regular)	9621	10
	Sterile Leg Bag 500ml variable inlet (regular)	9624	10
	Sterile Leg Bag 540ml with extension tubing and a pair of leg bag straps - latex free*	9348	12
	Sterile Leg Bag 800ml 10cm inlet (regular)	9631	10
	Sterile Leg Bag 800ml variable inlet (regular)	9632	10

*This product replaces the 9820 leg bag with a latex free option

Urinary Leg Bag Accessories		Order No.	Pack Size
	Leg Bag Straps	9625	10
	Deluxe Leg Bag Straps (35cm ex. small)	9663	2
	Deluxe Leg Bag Straps (38cm small)	9660	2
	Deluxe Leg Bag Straps (49cm medium)	9661	2
	Deluxe Leg Bag Straps (56cm large)	9662	2
	Leg Bag Holsters (small)	9612	4
	Leg Bag Holsters (medium)	9613	4
	Leg Bag Holsters (large)	9614	4
	Leg Bag Straps 14" 35.5cm (calf)	9342	1 pair
	Leg Bag Straps 23" 58.5cm (thigh) (for use with 9348 and 9814 Leg Bags)	9343	1 pair
	Universal Remover Wipes	7760	50
Skin Gel Protective Dressing Wipes	7917	50	

Drainage Bags 2000ml		Order No.	Pack Size
	Night Drainage Bags (Non-Drainable)	9651	10
	Night Drainage Bag (Drainable)	9650	10

For more information contact us FREE on:

UK 0800 521 377 Ireland 1800 503 400

For FREE samples contact us FREE on:

UK 0800 592 202 Email samples.uk@hollister.com
Ireland 1800 503 400 Email customerservices.ie@hollister.com

United Kingdom
Hollister Limited
Rectory Court
42 Broad Street
Wokingham
Berkshire RG40 1AB

Ireland
Hollister Limited
Unit 4045
Kingswood Road
Citywest Business Park
Co Dublin

www.hollister.co.uk